



Bed/Chair Occupancy Sensor

The bed/chair occupancy sensor will automatically trigger an alarm should it detect any of the following situations:

- If a client has not gone to bed by the predetermined time
- If a client has not got out of bed by the predetermined time
- If a client gets out of bed during the night and does not return within the predetermined time period

The solution can be tailored to suit the specific daily activities of the individual. There is an optional facility available to automatically turn on a bedside lamp or series of lights, such as those in the main bedroom, hallway and bathroom to minimize the risk of falling should the client get out of bed.

The bed/chair occupancy sensor can also be used as part of the virtual sensor feature on the Caresse+.



Features

- Wireless connection to the host unit (Caresse or Caresse+)
- Allows connection of a second bed pad for use in a double bed
- Automatic control of bedside lamp using optional X10 controller
- User can extend the absent time period using an optional “extend switch”
- Configurable for individual requirements using a PDA and Tunstall software

Technical Details

- Part Number:** 41004/13 (does not include the bed pad Tunstall part # D4106009A)
- Radio Frequency:** 312.000 MHz
- Battery:**
- > Four AA style Alkaline batteries (user replaceable)
 - > Typical life 12 months for 20 activations per day.
 - > Automatic low battery transmission to the host unit
- Dimensions:** 3.7"W x 6.2"H x 1.4" D (93 x 158 x 34mm)
- Weight:** 1.76 lbs (0.8Kg)
- Operating Temperature Range:** 14° F to 131° F (-10° C to +55° C)
- Radio Range:** 213 ft (65 m) in free space
- Supervision:** Full
- Standards:** FCC CFR47 part 15, Canadian RSS210
- Product Compatibility:** Caresse, Caresse+

Contact Details

Tunstall Group Limited
10N. Martingale Road Suite 400
PMB259
Schaumburg
Illinois 60173-2411
USA
Tel: 1-800-892-2205
Fax: +44 1977 660562
Email: sales@tunstallamerica.com

Tunstall Canada
7540 Bath Road
Mississauga
Ontario L4T 1L2
CANADA

Tel: 905-677-1144
Fax: 905-677-1121
Email: sales@tunstallamerica.com



www.tunstallamerica.com

Our policy of continual development means that product specification and appearance may change without notice. © 2008 Tunstall Group Ltd. ® TUNSTALL is a registered trade mark.