



Fall Detector

The Fall Detector is a unique personal device. It provides reliable detection of a potentially serious fall by the wearer and then triggers a call for help via a Tunstall Telecom personal emergency response unit. The device is attractively styled whilst being small and light enough to be worn unobtrusively.



Features

- Attractive and discreet
- Ergonomically styled with matte finish
- Manual trigger button
- Activation indicator light and audible sounder
- User replaceable battery
- Belt clip or holster wearing options
- Easy to keep clean
- 4 hour automatic transmission of status reports for pendant supervision
- Wireless connection to Telecom 400/4000 and Caresse/Caresse+
- Generates an audible warning as well as triggering a call on the host unit.

Technical Details

Part Number:	P66004/07 (light grey – as pictured) P66004/08 (charcoal grey)
Supplied parts:	Supplied in a box with user manual, battery and additional spare battery, belt and holster and table stand.
Radio Frequency:	312.000 MHz
Battery:	> User replaceable battery (Duracell® PX28L 6 volt) Typical life 4 to 5 months. > When the battery is nearing the end of its life, the light flashes when the button is pressed and the call recipient is automatically notified of the low battery state. Battery status is also reported automatically during the 4 hourly supervision transmissions
Dimensions:	3.0"W x 2.1"H x 1.1"D (75 x 53 x 28mm) Alarm button 1.1"W x 0.7"H (27x 17mm)
Weight:	0.17lbs (75g)
Operating Temperature Range:	14° F to 131° F (-10° C to +55° C)
Radio Range:	213 ft (65 m) in free space
Supervision:	Auto Low Battery
Standards:	FCC CFR47 part 15, Canadian RSS210
Product Compatibility:	Telecom 400, Telecom 4000, Caresse, Caresse+

Contact Details

Tunstall Group Limited
10N. Martingale Road Suite 400
PMB259
Schaumburg
Illinois 60173-2411
USA
Tel: 1-800-892-2205
Fax: +44 1977 660562
Email: sales@tunstallamerica.com

Tunstall Canada
7540 Bath Road
Mississauga
Ontario L4T 1L2
CANADA
Tel: 905-677-1144
Fax: 905-677-1121
Email: sales@tunstallamerica.com



www.tunstallamerica.com

Our policy of continual development means that product specification and appearance may change without notice. © 2008 Tunstall Group Ltd. ® TUNSTALL is a registered trade mark.