



## Gem Pendant

The Gem pendant is slightly larger in size than the Amie to provide a larger button and allow ease of use by people with even the most limited dexterity. Its attractive design and choice of fastening options make it a pleasure to wear.

It comes complete with a non-active neckcord and clothing clip which will clip to any type or thickness of material. It may also be worn with an optional wrist strap.

An additional easy press adaptor is also available to enable people with extremely limited dexterity to raise a call more easily.



## Features

- Attractive and discreet
- Ergonomically styled with smooth, gloss finish
- Activation light provides clear confirmation of button press and indication of low battery
- Water resistant design – IP67
- 4 hour automatic transmission of status reports for pendant supervision
- Range of easy wearing options (neck cord, wrist strap, clothing clip)
- Easy to keep clean

# Technical Details

- Part Number:** P64004/02 supplied in a box with neck cord and clothing clip
- Radio Frequency:** 312.000 MHz
- Battery:**
- > Integral battery (not user replaceable)
  - > Typical life 3.5 years for 20 activations for 20 activations per day, these include supervision transmissions
  - > When the battery is nearing the end of its life, the light flashes when the button is pressed and the call recipient is automatically notified of the low battery state. Battery status is also reported automatically during the 4 hourly supervision transmissions
- Dimensions:** 1.72" W x 1.87"H x 0.63"D (43.5 x 47.5 x 16mm)  
Button diameter 0.79" (20mm)
- Weight:** 0.06 lbs (25g)
- Operating Temperature Range:** 14° F to 131° F (-10° C to +55° C)
- Radio Range:** 295 ft (90 m) in free space
- Supervision:** Full
- Standards:** FCC CFR47 part 15, Canadian RSS210
- Product Compatibility:** Telecom 400 and Telecom 4000
- Optional Accessories:** Wrist strap part number D6402143A  
Easy Press Adaptor part number 64900/16A

# Contact Details

Tunstall Group Limited  
10N. Martingale Road Suite 400  
PMB259  
Schaumburg  
Illinois 60173-2411  
USA  
Tel: 1-800-892-2205  
Fax: +44 1977 660562  
Email: [sales@tunstallamerica.com](mailto:sales@tunstallamerica.com)

Tunstall Canada  
7540 Bath Road  
Mississauga  
Ontario L4T 1L2  
CANADA  
Tel: 905-677-1144  
Fax: 905-677-1121  
Email: [sales@tunstallamerica.com](mailto:sales@tunstallamerica.com)



[www.tunstallamerica.com](http://www.tunstallamerica.com)

Our policy of continual development means that product specification and appearance may change without notice. © 2008 Tunstall Group Ltd. ® TUNSTALL is a registered trade mark.